



106.3" Length | 34.5" Width | 30 Square Feet

AQUATIC THERAPY: Evolved to the size of your land-based treadmill

HydroWorx's most **ECONOMICAL** solution. This cutting-edge aquatic therapy system is designed to revolutionize rehabilitation and performance training. Featuring an underwater treadmill, resistance jets, and customizable water depths, it provides unmatched versatility for therapy, conditioning, and recovery. Engineered for durability and ease of use, the HydroWorx ECO enhances mobility, strength, and overall wellness in a safe, low-impact environment. Ideal for nearly any space with a footprint of only 30 square feet. The ECO maximizes patient outcomes and accelerates recovery at an affordable price.

The ultimate tool for providing top-notch care and achieving excellent therapeutic results.



Self-contained, no additional equipment required



Resistance jets & massage hose



Treadmill speed 0.2 - 6 mph



Fold-down seat



Designed for existing treatment spaces



Fits through a 36" door



Adjustable water depth




Build a profit with one 45-minute patient session/day

STANDARD FEATURES

- ◆ Self-contained, no additional equipment required
- ◆ Adjustable Water Depth from 31"-50"
- ◆ Treadmill Speed: 0.2-6 mph (in 0.1 mph increments)
- ◆ Resistance Jet with Deep Tissue Massage Hose
- ◆ Adjustable 7" Touch Screen User Controls
- ◆ Front and side handrails
- ◆ Fold Down Seat
- ◆ Chemical Feeder
- ◆ Adjustable Temperature Control

TECHNICAL SPECIFICATIONS

- ◆ Maximum Load: 300 lbs single user or 550 lbs with a Therapist
- ◆ Maximum water depth: 50"
- ◆ Fill/Empty Time: 1 min 30 seconds each way
- ◆ Temperature Range: ambient to 104 degrees
- ◆ Treadmill Tank plus Foyer Exterior Dimensions: 34.5" x 106.3"
- ◆ Complete System Weight:
 - Dry: 1500 lbs
 - With Water: 5200 lbs
- ◆ Capacity: 440 gallons
- ◆ Power Requirements: 240-volt, 40-amp, single phase



Medicare Reimbursement Model



One patient for a minimum 45 minute session (3 units)



97113 (Aquatic Therapy) pays approx. \$37 for 1st unit (\$27 after)



\$91/day x 21 business days in a month

= \$1,911 in Insurance Reimbursement Revenue
\$1,110 HydroWorx ECO monthly payment
and almost \$800 in profit